

Immunology Consultants Laboratory, Inc.

15862 SW 72nd Avenue, #150

Portland, OR 97224

T 503-747-2454

F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Purified Mouse IgG3 Kappa

CATALOG NUMBER: RS-90G3LY-K

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Lyophilized

BUFFER: PBS buffer with non-protein stabilizer.

STORAGE: Lyophilized material should be kept at 4°C until rehydrated and will remain stable for a minimum of 2 years. Rehydrate by adding 1.0 mL of reagent quality water to the product vial. Store the liquid solution at 4°C for up to a week; thereafter, aliquots should be prepared and stored frozen. Stable for a minimum of 1 year at -20°C.

PROCEDURE: IgG3 kappa hybridoma cells were grown in tissue culture media. The IgG3 was purified by a combination of salt fractionation and Protein A affinity chromatography.

SPECIFICATIONS: This product reacts only with anti-IgG3. It does not react with anti-IgG1, IgG2A, IgG2B, IgA, IgM and/or IgD by immunodiffusion and IEP techniques.

FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION

CERTIFICATE OF ANALYSIS VERIFICATION:
PRODUCT OF THE UNITED STATES OF AMERICA

DATE: _____

RS90G3LYK