

Immunology Consultants Laboratory, Inc.

15862 SW 72nd Avenue, #150

Portland, OR 97224

T 503-747-2454

F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Dog Fibrinogen Reference Standard

CATALOG NUMBER: RS-40FIB

LOT NUMBER:

QUANTITY: 1.0 ml

PHYSICAL STATE: Liquid

PRESERVATIVE: 0.1% Sodium Azide

STORAGE: 4C short term and frozen for >7 days.

PROCEDURE: Serum was collected from several dogs.

SPECIFICATIONS: Concentration and dilution determined by sandwich ELISA.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

CERTIFICATE OF ANALYSIS VERIFICATION:
PRODUCT OF THE UNITED STATES OF AMERICA

DATE: