



## **PRODUCT PROFILE**

**PRODUCT:** Affinity Purified Goat anti-Mouse C3

**CATALOG NUMBER:** GC3-90A-Z

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Goats were immunized with highly purified mouse C3 and the resulting antiserum was collected. Antibody was immunoaffinity purified off an C3 containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with mouse C3 as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with C3 from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***