

PRODUCT PROFILE

PRODUCT: Purified Hamster (CHO) Clusterin

CATALOG NUMBER: AG65-0075-01-Z

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS

PRESERVATIVE: 0.09% sodium azide

STORAGE: Store at -20°C. The thawed product is stable at 2-8°C for one week. Following the initial thaw, please prepare multiple aliquots and store at -20°C. It is not recommended to store the product in a frost-free freezer. Multiple freeze-thaw cycles are not recommended.

EXPIRATION: 12 months from shipment date.

PROCEDURE: Recombinant CHO Clusterin (GenBank EGW05785.1) with a 6xHIS tag was expressed in CHO-S cells and purified by nickel affinity chromatography.

SPECIFICATIONS: The product is determined to be >90% pure based on SDS-PAGE analysis. Under non-reducing conditions, clusterin is detected at an estimated molecular weight of 60-70kDa. Its concentration is determined using a Bradford protein assay. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY. NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY. NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

PRODUCT OF THE UNITED STATES OF AMERICA

AG65-0075-01-Z_Rev021924