

PRODUCT PROFILE

PRODUCT: Affinity Purified Chicken anti-*Campylobacter jejuni*

CATALOG NUMBER: CCBJ-55A-202

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Chickens were immunized with highly purified *Campylobacter jejuni* hippurate hydrolase (HipO) and the eggs were collected and delipidated. Antibody was immunoaffinity purified off a HipO containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody reacts with *Campylobacter jejuni* hippurate hydrolase as determined by ELISA and it can detect *Campylobacter jejuni* cells with Immunofluorescence staining. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA