

## **PRODUCT PROFILE**

**PRODUCT:** Affinity Purified Chicken anti-Canine Leptospira LipL32 Antibody

**CATALOG NUMBER:** CL32-40A

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Chickens were immunized with highly purified canine Leptospira LipL32 protein and the eggs were collected and delipidated. Antibody was immunoaffinity purified off a canine Leptospira LipL32 protein containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

**SPECIFICATIONS:** This antibody reacts with canine Leptospira LipL32 protein and is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***