

PRODUCT PROFILE

PRODUCT: Affinity Purified Chicken anti-Marburg Virus VP40

CATALOG NUMBER: CM40-55A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Chickens were immunized with highly purified Marburg Virus VP40 and the eggs were collected and delipidated. Antibody was immunoaffinity purified off a Marburg Virus VP40 containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody reacts with Marburg Virus VP40 as determined by ELISA techniques. It is suitable for blotting and ELISA applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA