

## **PRODUCT PROFILE**

- PRODUCT:** Affinity Purified Goat anti-Cholera Toxin A+B
- CATALOG NUMBER:** GCHL-55A-Z
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** PBS, pH 7.2
- PRESERVATIVE:** 0.09% sodium azide
- STORAGE:** 2 - 8° C
- EXPIRATION:** 24 months from shipping date.
- PROCEDURE:** Goats were immunized with highly purified intact vibrio cholera toxin and the resulting antiserum was collected. Antibody was immunoaffinity purified off a cholera toxin containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
- SPECIFICATIONS:** This antibody reacts with intact cholera toxin, the alpha subunit and the beta subunit as determined by ELISA. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

GCHL-55A-Z\_Rev121223