

## **PRODUCT PROFILE**

**PRODUCT:** Purified Mouse anti-BVDV (Clone #104)

**CATALOG NUMBER:** MBVD-10A-104

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** Antibody should be stored between 2°C and 8°C for up to 30 days . For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.

**EXPIRATION:** 12 months from shipping date.

**PROCEDURE:** This mouse IgG1 antibody is purified using salt fractionation procedures and Protein A chromatography. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with Bovine Viral Diarrhea Virus as determined by ELISA and IFA. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY; NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

***PRODUCT OF THE UNITED STATES OF AMERICA***

*MBVD-10A-104\_Rev042624*