

PRODUCT PROFILE

- PRODUCT:** Purified Mouse anti-FeLV (Clone #CM1)
- CATALOG NUMBER:** MFLV-20M-CM1
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** PBS, pH 7.2
- PRESERVATIVE:** 0.09% Sodium Azide
- STORAGE:** Store between 2 and 8 °C for up to 30 days. For long-term storage, aliquot and store between -20 °C and -80 °C. Multiple freeze-thaw cycles are not recommended.
- EXPIRATION:** 12 months from shipping date.
- PROCEDURE:** This mouse IgG1 is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
- SPECIFICATIONS:** This antibody reacts with feline leukemia virus (FeLV) p27 as determined by ELISA, LFD and western blot. Suggest using this as ELISA capture antibody when detecting with MFLV-20M-8D11. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR CAT OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA