

PRODUCT PROFILE

- PRODUCT:** Purified Mouse anti-Canine Procalcitonin (Clone #11F9)
- CATALOG NUMBER:** MPCT-40A-11F9
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** PBS, pH 7.2
- PRESERVATIVE:** 0.09% Sodium Azide
- STORAGE:** Antibody should be stored between 2 and 8°C for up to 14 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
- EXPIRATION:** 12 months from shipping date.
- PROCEDURE:** This mouse IgG2a kappa is purified using Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
- SPECIFICATIONS:** This antibody reacts with canine procalcitonin as determined by ELISA and western blotting. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA