

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT:	Purified Annexin A5 CHO-S
CATALOG NUMBER:	AG65-0425-01-Z
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS
PRESERVATIVE:	0.09% sodium azide
STORAGE:	Store at -20°C. The thawed product is stable at 2-8°C for one week. Following the initial thaw, please prepare multiple aliquots and store at -20°C. It is not recommended to store the product in a frost-free freezer. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipment date.
PROCEDURE:	Recombinant CHO Annexin A5 (GenBank XP_027247784.1) with a 6xHIS tag was expressed in CHO-S cells and purified by nickel affinity chromatography.
SPECIFICATIONS:	The product is determined to be >90% pure based on SDS-PAGE analysis. Under reducing conditions, Annexin A5 is detected at an estimated molecular weight of 40kDa. Its concentration is determined using a Bradford protein assay. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY. NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY. NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

PRODUCT OF THE UNITED STATES OF AMERICA