

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Affinity Purified Chicken anti-Green Fluorescent Protein

CATALOG NUMBER: CGFP-45A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Chickens were immunized with highly purified green-fluorescent protein (GFP) isolated directly from *Aequorea victoria*. Eggs were collected. Antibodies were immunoaffinity purified using the native GFP immobilized on a solid phase.

SPECIFICATIONS: ICL's anti-GFP recognizes some variants of recombinant *Aequorea* green-fluorescent protein. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

CGFP-45A_Rev031622