



PRODUCT PROFILE

PRODUCT: Affinity Purified Chicken anti-Rabbit IgG h+l – cross adsorbed antibody
Minimum reactivity with human and mouse.

CATALOG NUMBER: CGHL-15AX8090

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Chickens were immunized with highly purified rabbit IgG h+l and eggs were collected. Antibody was solid phase adsorbed using human and mouse immunosorbents to remove cross reactive antibody. Specific antibody was immunoaffinity purified off a rabbit IgG containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody will react with rabbit IgG and with light chains common to other rabbit immunoglobulins. It has less than 1% cross reactivity with human and mouse IgG as determined by ELISA and IEP techniques. This antibody is suitable for blotting, ELISA and IHC applications. This antibody may cross react with IgG from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA