

## Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

### **PRODUCT PROFILE**

<b>PRODUCT:</b>	<b>Purified Mouse anti-Coronavirus (SARS-CoV-2) - Clone# 111</b>
<b>CATALOG NUMBER:</b>	MCV19-80A-111
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	1.0 mg
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, pH 7.2.
<b>PRESERVATIVE:</b>	0.09% sodium azide
<b>STORAGE:</b>	Antibody should be stored between 2 – 8 °C for up to 180 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
<b>EXPIRATION:</b>	12 months from shipping date.
<b>PROCEDURE:</b>	This mouse IgG1 is purified using salt fractionation procedures and Protein A chromatography.
<b>SPECIFICATIONS:</b>	This antibody reacts with the nucleoprotein of SARS-CoV-2 as determined by ELISA, western blot and lateral flow. Optimal working dilutions should be determined experimentally by the investigator. There may be some cross-reactivity with other strains of coronavirus.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

MCV19-80A-111\_Rev102020