

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223

T 503-747-2454

F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Purified Mouse anti-Human Ferritin (Clone #101)

CATALOG NUMBER: MF-80A-101

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2.

PRESERVATIVE: 0.1% sodium azide

STORAGE: Antibody should be stored between 2 – 8 °C for up to 180 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.

EXPIRATION: 12 months from shipping date.

PROCEDURE: This mouse IgG1 is purified using Protein A chromatography.

SPECIFICATIONS: This antibody reacts with human ferritin as determined by ELISA and lateral flow. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

MF-80A-101_Rev072121