



PRODUCT PROFILE

PRODUCT:	Purified Mouse anti-Bovine Haptoglobin (Clone #1B5)
CATALOG NUMBER:	MHPT-10A-1B5
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.09% sodium azide
STORAGE:	Antibody should be stored between 2°C and 8°C for up to 14 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipping date.
PROCEDURE:	This mouse IgG1 kappa antibody is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with recombinant E2 protein from classical swine fever virus (CSFV) as determined by ELISA, LFD and western blot. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY; NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

PRODUCT OF THE UNITED STATES OF AMERICA