



## **PRODUCT PROFILE**

<b>PRODUCT:</b>	Purified Mouse anti- <i>Listeria monocytogenes</i> (Clone #6D5)
<b>CATALOG NUMBER:</b>	MLSTA-55A-6D5
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, pH 7.2
<b>PRESERVATIVE:</b>	0.09% sodium azide
<b>STORAGE:</b>	Antibody should be stored between 2°C and 8°C for up to 14 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
<b>EXPIRATION:</b>	12 months from shipping date.
<b>PROCEDURE:</b>	This mouse IgG1 kappa antibody is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
<b>SPECIFICATIONS:</b>	This antibody reacts with <i>Listeria monocytogenes</i> as determined by ELISA and western blotting. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY; NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

**PRODUCT OF THE UNITED STATES OF AMERICA**