



PRODUCT PROFILE

PRODUCT:	Purified Mouse anti-Human Procalcitonin (Clone #10F3)
CATALOG NUMBER:	MPCT-80A-10F3
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	Phosphate buffered saline, pH 7.2
PRESERVATIVE:	0.09% Sodium Azide
STORAGE:	Antibody should be stored between 2 and 8°C for up to 14 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipping date.
PROCEDURE:	This mouse IgG1 kappa is purified using Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with Human Procalcitonin (hPCT) as determined by ELISA and western blotting. Its epitope specificity has been mapped to the Calcitonin fragment of hPCT. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA