



PRODUCT PROFILE

PRODUCT:	Purified Mouse anti-Prostate Specific Antigen (PSA - Clone #803)
CATALOG NUMBER:	MPSA-80A-803
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2.
PRESERVATIVE:	0.09% sodium azide
STORAGE:	Antibody should be stored between 2 – 8 °C for up to 180 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipping date.
PROCEDURE:	This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with prostate specific antigen as determined by ELISA and western blotting. Optimal working dilutions should be determined experimentally by the investigator. This clone has been demonstrated to pair with clone #801.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY; NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

PRODUCT OF THE UNITED STATES OF AMERICA