



PRODUCT PROFILE

PRODUCT:	Purified Mouse anti-Phospho Tau-217 (Clone #100)
CATALOG NUMBER:	MT217-80A-100
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.09% sodium azide
STORAGE:	Antibody should be stored between 2°C and 8°C for up to 30 days . For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipping date.
PROCEDURE:	This mouse IgG1 antibody is purified using salt fractionation procedures and Protein A chromatography. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.
SPECIFICATIONS:	This antibody reacts with pTau217 as determined by ELISA and LFD. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY; NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

PRODUCT OF THE UNITED STATES OF AMERICA