



PRODUCT PROFILE

PRODUCT:	Purified Mouse anti-Thyroxine (T4) (Clone #101)
CATALOG NUMBER:	MT4-65A-101
LOT NUMBER:	
QUANTITY:	
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.09% sodium azide
STORAGE:	Antibody should be stored between 2°C and 8°C for up to 30 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
EXPIRATION:	12 months from shipping date.
PROCEDURE:	This mouse IgG1 antibody is purified using salt fractionation procedures and Protein A chromatography.
SPECIFICATIONS:	This antibody reacts with Thyroxine (T4) as determined by ELISA and LFD. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY; NOT FOR HUMAN OR ANIMAL CONSUMPTION.**

PRODUCT OF THE UNITED STATES OF AMERICA