

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544  
[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## PRODUCT PROFILE

**PRODUCT:** Mouse IgG Enriched

**CATALOG NUMBER:** RA-90G

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% Sodium Azide

**STORAGE:** 2 – 8°C

**EXPIRATION:** 12 months from shipment date

**PROCEDURE:** The Ig fraction was enriched from normal mouse serum using a combination of procedures that include salt fractionation.

**SPECIFICATIONS:** Antibody concentration was determined by extinction coefficient: O.D. 1.4 at 280nm equals 1.0 mg of IgG. Estimated purity of >80% as determined by SDS-PAGE. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

RA-90G\_Rev043020