

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

- PRODUCT:** Rabbit anti-BSA - Antiserum
- CATALOG NUMBER:** RAL-10
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- PRESERVATIVE:** 0.09% sodium azide
- STORAGE:** 2 - 8° C
- EXPIRATION:** 24 months from shipping date.
- PROCEDURE:** Rabbits were immunized with highly purified bovine serum albumin (BSA) and the resulting antiserum was collected. The antiserum was delipidated, solid phase adsorbed and 0.45 um filtered. Total protein measurement was determined by refractometer.
- SPECIFICATIONS:** This antiserum reacts with BSA as determined by IEP and Ouchterlony techniques. This antiserum may cross react with albumin from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR BOVINE OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA