



PRODUCT PROFILE

PRODUCT: Biotin Conjugated Rabbit anti-CBP Tag

CATALOG NUMBER: RCBP-45B

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2 with stabilizer

PRESERVATIVE: 0.09% Sodium Azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Rabbits were immunized with highly purified MKRRWKKNFIAVSAANRFKKISSSGAL, Calmodulin Binding Protein (CBP) tag, conjugated with KLH and the resulting antiserum was collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent and conjugated with biotin. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with MKRRWKKNFIAVSAANRFKKISSSGAL as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA