

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544  
[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## PRODUCT PROFILE

**PRODUCT:** Normal Rabbit Serum – Sterile Filtered

**CATALOG NUMBER:** RS-15-100ML

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** None Added

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** Material as received can be kept at 2 - 8° C for up to 10 days. Aliquot and store between -20°C and -80°C for greater than 10 days.

**PROCEDURE:** Normal rabbit serum was collected from healthy rabbits and sterile filtered.

**SPECIFICATIONS:** Optimal working dilutions should be determined experimentally by the investigator.

**Total Protein:** \_\_\_\_\_ mg/mL

FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION

*PRODUCT OF THE UNITED STATES OF AMERICA*