

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Normal Cat Serum – Sterile Filtered

CATALOG NUMBER: RS-20-100ML

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: None Added

PRESERVATIVE: 0.09% sodium azide

STORAGE: Material as received can be kept at 2 - 8° C for up to 10 days. Aliquot and store between -20°C and -80°C for greater than 10 days.

PROCEDURE: Normal cat serum was collected from healthy cats and sterile filtered.

SPECIFICATIONS: Optimal working dilutions should be determined experimentally by the investigator.

Total Protein: _____ mg/mL

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

RS-20-100ML_Rev030222