

## Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

### PRODUCT PROFILE

**PRODUCT:** Rat IgG Purified

**CATALOG NUMBER:** AG-25G

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipment date

**PROCEDURE:** The IgG h+I was purified from rat serum using a combination of procedures that include salt fractionation and protein A purification.

**SPECIFICATIONS:** Antibody concentration was determined by extinction coefficient: O.D. 1.40 at 280nm equals 1.0 mg of IgG. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

AG-25G\_Rev050120