

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223

T 503-747-2454

F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Dog IgG Purified

CATALOG NUMBER: AG-40G

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipment date

PROCEDURE: The IgG h+I was purified from dog serum using a combination of procedures that include salt fractionation and protein A purification.

SPECIFICATIONS: Antibody concentration was determined by extinction coefficient: O.D. 1.40 at 280nm equals 1.0 mg of IgG. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

AG-40G_Rev050120