



PRODUCT PROFILE

PRODUCT: HRP Conjugated Chicken anti-Cotton Rat IgG h+l

CATALOG NUMBER: CCOT-25P

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: HEPES – saline with stabilizer

PRESERVATIVE: 0.002% thimerosal

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Chickens were immunized with highly purified Cotton rat IgG and eggs were collected. Antibody was sold phase adsorbed to remove cross reactive antibody. Specific antibody was immunoaffinity purified off a Cotton Rat IgG containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody reacts with Cotton Rat IgG as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with IgG from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA