



## **PRODUCT PROFILE**

<b>PRODUCT:</b>	Affinity Purified Chicken anti-Digoxigenin
<b>CATALOG NUMBER:</b>	CDIG-65A
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, pH 7.2
<b>PRESERVATIVE:</b>	0.09% sodium azide
<b>STORAGE:</b>	2 - 8° C
<b>EXPIRATION:</b>	24 months from shipment date
<b>PROCEDURE:</b>	Chickens were immunized with digoxigenin and egg yolks were collected. An IgY enriched fraction was prepared, the antibody was immunoaffinity purified off an antigen containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.
<b>SPECIFICATIONS:</b>	This antibody reacts with digoxigenin as determined by blotting and ELISA techniques. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***