



PRODUCT PROFILE

| | |
|------------------------|--|
| PRODUCT: | Affinity Purified Chicken anti-DYKDDDDK |
| CATALOG NUMBER: | CFLG-45A-Z |
| LOT NUMBER: | |
| QUANTITY: | |
| PHYSICAL STATE: | Liquid |
| BUFFER: | PBS, pH 7.2 |
| PRESERVATIVE: | 0.09% sodium azide |
| STORAGE: | 2 - 8° C |
| EXPIRATION: | 24 months after shipping date. |
| PROCEDURE: | Chickens were immunized with the DYKDDDDK peptide conjugated to KLH. Eggs were collected and the IgY fraction prepared. Antibody was immunoaffinity purified using peptide coupled to a different carrier. Antibody concentration was determined using an absorbance at 290 nm: 1.32 equals 1.0 mg of IgY. |
| SPECIFICATIONS: | Optimal working dilutions should be determined for your particular assay conditions. |

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA