



PRODUCT PROFILE

PRODUCT: Affinity Purified Chicken anti beta-Galactosidase

CATALOG NUMBER: CGAL-45A-Z

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Chickens were immunized with highly purified beta-galactosidase from *Escherichia coli*. Eggs were collected and an IgY fraction prepared. Antibodies were solid phase absorbed then immunoaffinity purified using purified beta-galactosidase immobilized on a solid phase. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody reacts with beta-galactosidase as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA