



## **PRODUCT PROFILE**

**PRODUCT:** Affinity Purified Chicken anti-HA Tag

**CATALOG NUMBER:** CHGT-45A

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Chickens were immunized with YPYDVPDYA (influenza hemagglutinin-HA-epitope). Eggs were collected and an IgY enriched fraction was prepared. Antibody was solid phase adsorbed to ensure its specificity. The anti-HA tag antibody was immunoaffinity purified off an antigen containing immunoabsorbant.

**SPECIFICATIONS:** This antibody reacts with YPYDVPDYA as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***