

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: HRP Conjugated Chicken anti-6xHIS Tag

CATALOG NUMBER: CHIS-45P-Z

LOT NUMBER:

QUANTITY:

BUFFER: 50 mM HEPES – Saline with Stabilizer

PRESERVATIVE: 0.002% Thimerosal

STORAGE: 2 - 8° C

EXPIRATION: 12 months after shipping date.

PROCEDURE: Chickens were immunized with highly purified HHHHHH conjugated with KLH (6xHIS). Eggs were collected. Antibody was immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody reacts with HHHHHH as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

CHIS-45P-Z_Rev082020