



PRODUCT PROFILE

PRODUCT: Biotin Conjugated Chicken anti-C-MYC

CATALOG NUMBER: CMYC-45B

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2 with Stabilizer

PRESERVATIVE: 0.09% Sodium Azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipment date

PROCEDURE: Chickens were immunized with EQKLISEEDL (C-MYC) conjugated with KLH. Eggs were collected and an IgY enriched fraction was prepared. Antibody was solid phase adsorbed to ensure its specificity before it was immunoaffinity purified off an antigen containing immunoabsorbant and conjugated with biotin. Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

SPECIFICATIONS: This antibody reacts with EQKLISEEDL as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA