

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## **PRODUCT PROFILE**

**PRODUCT:** HRP Conjugated Chicken anti-C-MYC

**CATALOG NUMBER:** CMYC-45P

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** 50 mM HEPES – Saline with Stabilizer

**PRESERVATIVE:** 0.002% Thimerosal

**STORAGE:** 4° C

**EXPIRATION:** 12 months after shipping date.

**PROCEDURE:** Chickens were immunized with EQKLISEEDL (C-MYC). Eggs were collected and IgY fraction prepared. Antibody was solid phase adsorbed, immunoaffinity purified off a peptide containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.

**SPECIFICATIONS:** This antibody reacts with EQKLISEEDL as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

CMYC-45P\_Rev041019