



## **PRODUCT PROFILE**

- PRODUCT:** HRP Conjugated Chicken anti-Protein A
- CATALOG NUMBER:** CPA-65P
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** HEPES – saline with stabilizer
- PRESERVATIVE:** 0.002% thimerosal
- STORAGE:** 2 - 8° C
- EXPIRATION:** 12 months from shipping date.
- PROCEDURE:** Chickens were immunized with recombinant Protein A and egg yolks were collected. An IgY enriched fraction was prepared, antibody was immunoaffinity purified off an antigen containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.
- SPECIFICATIONS:** This antibody reacts with Protein A as determined by blotting and ELISA techniques. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***