

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## PRODUCT PROFILE

- PRODUCT:** HRP Conjugated Chicken anti-Protein G
- CATALOG NUMBER:** CPG-65P
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** HEPES – saline with stabilizer
- PRESERVATIVE:** 0.002% thimerosal
- STORAGE:** 2 - 8° C
- EXPIRATION:** 12 months from shipping date.
- PROCEDURE:** Chickens were immunized with highly purified Protein G and the resulting antiserum was collected. Antibody was solid phase adsorbed, immunoaffinity purified off a Protein G containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgG.
- SPECIFICATIONS:** This antibody reacts with Protein G as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

CPG-65P\_Rev101819