



## **PRODUCT PROFILE**

- PRODUCT:** Affinity Purified Chicken anti-Protein L
- CATALOG NUMBER:** CPL-65ALY
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Lyophilized
- BUFFER:** PBS with non-protein stabilizer
- PRESERVATIVE:** None added
- STORAGE:** Lyophilized material should be kept at 4°C until rehydrated and will remain stable for a minimum of 2 years.  
Rehydrate by adding 1.0 mL of reagent quality water to the product vial and add preservative. Liquid solution stored at 2-8°C is stable for a minimum of 1 year.
- EXPIRATION:** 24 months from shipment date (lyophilized form).
- PROCEDURE:** Chickens were immunized with recombinant Protein A and egg yolks were collected. An IgY enriched fraction was prepared, the antibody was immunoaffinity purified off an antigen containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.32 equals 1.0 mg of IgY.
- SPECIFICATIONS:** This antibody reacts with Protein A as determined by blotting and ELISA techniques. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

**PRODUCT OF THE UNITED STATES OF AMERICA**