



PRODUCT PROFILE

PRODUCT: Affinity Purified Goat anti-Chicken Avian Leukosis Virus p27

CATALOG NUMBER: GALV-30A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Goats were immunized with purified recombinant avian leukosis virus (ALV) p27 and the resulting antiserum was collected. Antibody was immunoaffinity purified off an ALV p27 containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with recombinant ALV p27 as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA, lateral flow and IHC applications. This antibody may cross react with p27 protein from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR BOVINE OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA