

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: HRP Conjugated Goat anti-Biotin

CATALOG NUMBER: GBIO-65P

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2 with stabilizer

PRESERVATIVE: 0.002% Thimerosal

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Goats were immunized with Biotin-KLH and the resulting antiserum was collected. Antibody was immunoaffinity purified off an antigen containing immunosorbent and conjugated with horseradish peroxidase (HRP). Antibody concentration was determined prior to conjugation using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with biotin as determined by blotting and ELISA techniques. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA

GBIO-65P_Rev030819