



PRODUCT PROFILE

PRODUCT: FITC Conjugated Goat anti-Mouse C3

CATALOG NUMBER: GC3-90F-Z

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2 with 0.2% BSA.

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 12 months from shipping date.

PROCEDURE: Goats were immunized with highly purified mouse C3 and the resulting antiserum was collected. Antibody was immunoaffinity purified off a mouse C3 containing immunosorbent and conjugated with fluorescein isothiocyanate (FITC). Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG. FITC/Protein Ratio = 3.0 – 6.0.

SPECIFICATIONS: This antibody will react with mouse C3 as determined by ELISA and IEP techniques. It is suitable for blotting, ELISA and IF applications. This antibody may cross react with C3 from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA