

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223

T 503-747-2454

F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Affinity Purified Goat anti-Human C4

CATALOG NUMBER: GC4-80A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Goats were immunized with highly purified human C4 and the resulting antiserum was collected. Antibody was immunoaffinity purified off a C4 containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with human C4 as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. This antibody may cross react with C4 from other species. Optimal working dilutions should be determined experimentally by the investigator.

FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION

PRODUCT OF THE UNITED STATES OF AMERICA

GC4-80A_Rev030119