



PRODUCT PROFILE

PRODUCT: Goat anti-Dog CRP - Antiserum

CATALOG NUMBER: GCRP-40

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Goats were immunized with highly purified dog c-reactive protein (CRP) and the resulting antiserum was collected. The antiserum was delipidated, solid phase adsorbed and 0.45 um filtered. Total protein measurement was determined by refractometer. Specific antibody content is determined immunoaffinity purification off an antigen containing immunosorbent.

SPECIFICATIONS: This antiserum reacts with dog CRP as determined by IEP and Ouchterlony techniques. This antiserum may cross react with CRP from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA