

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Affinity Purified Goat anti-Horse Hemoglobin

CATALOG NUMBER: GHM-70A

LOT NUMBER:

QUANTITY:

PHYSICAL STATE: Liquid

BUFFER: PBS, pH 7.2

PRESERVATIVE: 0.09% sodium azide

STORAGE: 2 - 8° C

EXPIRATION: 24 months from shipping date.

PROCEDURE: Goats were immunized with highly purified horse hemoglobin and the resulting antiserum was collected. Antibody was immunoaffinity purified off a hemoglobin containing immunosorbent. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

SPECIFICATIONS: This antibody reacts with horse hemoglobin as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA, lateral flow and IHC applications. This antibody may cross react with hemoglobin from other species. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HORSE OR ANIMAL CONSUMPTION**

PRODUCT OF THE UNITED STATES OF AMERICA