

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## **PRODUCT PROFILE**

**PRODUCT:** Goat anti-Rabbit Whole Serum

**CATALOG NUMBER:** GWS-15

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** N/A

**PRESERVATIVE:** 0.09% sodium azide

**STORAGE:** 2 - 8° C

**EXPIRATION:** 24 months from shipping date.

**PROCEDURE:** Goats were immunized with normal rabbit serum and the resulting antiserum was collected. Antibody concentration was determined using a refractometer.

**SPECIFICATIONS:** This antibody reacts with multiple proteins in rabbit serum when tested by immunoelectrophoresis (IEP). Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR RABBIT OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

GWS-15\_Rev031419