

Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street
Portland, OR 97223
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

- PRODUCT:** Purified Mouse anti-Alexa Fluor 647
- CATALOG NUMBER:** M647-65A-400
- LOT NUMBER:**
- QUANTITY:**
- PHYSICAL STATE:** Liquid
- BUFFER:** PBS, pH 7.2.
- PRESERVATIVE:** 0.09% sodium azide
- STORAGE:** Antibody should be stored between 2 – 8 °C for up to 14 days. For long-term storage, aliquot and store between -20°C and -80°C. Multiple freeze-thaw cycles are not recommended.
- EXPIRATION:** 12 months from shipping date.
- PROCEDURE:** This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography.
- SPECIFICATIONS:** This antibody reacts with Alexa Fluor 647 conjugated proteins as determined by ELISA, western blot and lateral flow. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

M647-65A_Rev032720