

# Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223  
T 503-747-2454  
F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

## **PRODUCT PROFILE**

**PRODUCT:** Purified Mouse anti-Cat IgE (Clone #2C1)

**CATALOG NUMBER:** ME-20A-2C1

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** PBS, pH 7.2

**PRESERVATIVE:** 0.09% Sodium Azide

**STORAGE:** 2 - 8° C for less than 30 days. Aliquot and store between -20C to -80C for long term storage. Avoid multiple freeze-thaw cycles.

**EXPIRATION:** 12 months from shipping date.

**PROCEDURE:** This mouse IgG2b kappa is purified using salt fractionation procedures and Protein A chromatography. It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with cat IgE (epsilon chain specific) as determined by ELISA and western blot. It is suitable for use in blotting assays, ELISA and lateral flow. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR CAT OR ANIMAL CONSUMPTION**

***PRODUCT OF THE UNITED STATES OF AMERICA***

ME-20A-2C1\_Rev031319