

## Immunology Consultants Laboratory, Inc.

7150 SW Sandburg Street  
Portland, OR 97223

T 503-747-2454

F 503-747-2544

[www.icllab.com](http://www.icllab.com)



ANTIBODIES ARE OUR BUSINESS

### **PRODUCT PROFILE**

**PRODUCT:** HRP Conjugated Mouse anti-Human IgE (Clone 24A)

**CATALOG NUMBER:** ME-80P-24A

**LOT NUMBER:**

**QUANTITY:**

**PHYSICAL STATE:** Liquid

**BUFFER:** HEPES - saline with stabilizer

**PRESERVATIVE:** 0.002% thimerosal

**STORAGE:** 2 - 8° C

**EXPIRATION:** 12 months from shipping date.

**PROCEDURE:** This mouse IgG1 kappa is purified using salt fractionation procedures and Protein A chromatography and conjugated with horseradish peroxidase (HRP). It is produced from tissue culture media obtained at our facility in The United States. Antibody concentration was determined using an absorbance at 280 nm: 1.4 equals 1.0 mg of IgG.

**SPECIFICATIONS:** This antibody reacts with human IgE (epsilon chain specific) as determined by ELISA and western blot. It is suitable for use in blotting assays and ELISA. Optimal working dilutions should be determined experimentally by the investigator.

**FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.  
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION**

**PRODUCT OF THE UNITED STATES OF AMERICA**